FOBT


The FOBT (Faecal occult blood test) tab provides three subtabs:

FOBT Last Recorded

A pie chart displays the % of patients with a FOBT Last Recorded in the following categories: 0-2 years, >2-3 years, >3-4 years, > 4 years, No FOBT Recorded, Ineligible*.

(1)

This report needs filtering for the appropriate age group as by default it will show all patients!




• FOBT Orders and results


A pie chart displays the % of patients who are required to have an FOBT test (i.e. are not ineligible) and have had an FOBT test order or result in the last 12 months in the following categories:

- FOBT Order and result an FOBT test result has been received in the last 12 months with an order recorded within 6 months prior
- FOBT Order only an FOBT order has been recorded but no result has been received
- FOBT Result only an FOBT order has been received but no order is recorded

This graph allows the clinician to follow up patients who have had an order placed but where a result has not been received.



The third tab, FOBT Summary Report Card, shows a summary of the available data broken down by all patients and active patients by age groups and time of test done.



FOBT Mappings with clinical systems

The FOBT mappings are provided in the 'General Data Mappings' Appendices for each clinical application.

*Patients with certain coded diagnosis are excluded from this report and listed under ineligible:

Excluded from FOBT

Anal carcinoma

Cancer of the anus

Carcinoma of the anus

Squamous cell carcinoma of anus

Adenocarcinoma of the colon

Bowel cancer

Cancer of the colon

Carcinoma of the colon

Colon carcinoma

Colonic cancer

Colonic carcinoma

Cancer of the rectum

Carcinoma of the rectum

Rectal adenocarcinoma

Rectal cancer

Rectal carcinoma

Genetic tendency to Hereditary nonpolyposis colorectal cancer

Hereditary nonpolyposis colorectal carcinoma

Lynch syndrome

Carcinosarcoma of the colon

Colon carcinosarcoma

Signet ring cell carcinoma of bowel


Signet ring cell carcinoma of rectum

Mammogram/FOBT for MD Users


Mammogram and FOBT results are not sent electronically in all states, but if test names are entered manually into the results tab CAT4 will pick up that the test has been done. The following test names are recognised:


Mammogram test names	FOBT test names
Breast Mammogram Screening Bilateral Mammography Ultrasound Breast Bilateral Wesley Breast Clinic Consultation Report Mammogram Mammogram-normal Breast Mammography Mammography - Bilateral Mammogram Breastscreen NSW Screening Mammogram - BreastScreen Queensland	IFOBT FOBT many others as in the 'general data mapping' section of your CIS

To add a test name manually to the Results tab in MD and right click in the results area to bring up the 'Add' option:



Then add one of the test names above and save the result (there has to be some text in the result field but CAT4 will only extract the test name from the 'Subject' field:





MD: General Data Category Mappings MD


Mammogram/FOBT for BP users

Mammogram and FOBT results are not sent electronically in all states, but if test names are entered manually into the results tab CAT4 will pick up that the test has been done. The following test names are recognised:

Mammogram test names	FOBT test names
----------------------	-----------------

Breast Mammogram Screening
Bilateral Mammography
Ultrasound Breast Bilateral
Wesley Breast Clinic Consultation Report
Mammogram
Mammogram-normal
Breast Mammography
Mammography - Bilateral
Mammogram Breastscreen NSW
Screening Mammogram - BreastScreen Queensland

To manually add a mammogram result click on the 'Add' button under the 'Investigation report' tab, then enter one of the test names from above and complete the other fields. CAT4 will only extract the test name and date but no outcome/result.



BP: General Data Category Mappings BP