

# After Server Installation

- By default, the service will be installed to:
  - **32Bit** – 'C:\Program Files\PCS\PCS Clinical Service'
  - **64Bit** – 'C:\Program Files (x86)\PCS\PCS Clinical Service'
- Start Menu shortcuts will be registered in 'PCS ->PCS Clinical Service'
  - **'Stop Service'** – Stops the PEN CS Clinical Service
  - **'Start Service'** – Starts the PEN CS Clinical Service
  - **'Runtime log'** – Short cut to the log file (same as Sidebar's)
  - **'Installation log'** – Shortcut to the installation log
- Add/Remove programs entry added 'PCS Clinical Service' (Can be un-installed from here)
- "PCS Clinic Service" and "PCS Clinic Monitoring Service" services should be running
- Firewall exception should be added to port 8086 (default) so that client installs can talk to the server.



Best Practice users please note: If you are running Topbar and you are trying to restore a BP Backup, the running Topbar service will prevent the backup from being restored. Please stop the PCS Clinical Service and PCS Clinical Monitoring Service as described above before attempting a database restore, then restart the service!

To install the client module on individual workstations please follow the steps in the [Client Module installation](#) guide