

Client installation steps

Go to <https://topbar-users.pencs.com.au>

Log on using the username and password provided by PENCS or created by your clinic administrator


The image shows a 'Sign In' form. At the top, the text 'Sign In' is centered. Below it are two input fields: the first is labeled 'User Name' with a person icon, and the second is labeled 'Password' with a lock icon. A blue 'LOGIN' button is positioned below the password field. At the bottom left, there is a link that says 'Forgotten your password?' with a question mark icon.

Image 1: Topbar logon screen

Click on the clinic you are installing to




Image 2: Topbar clinic screen

Click on the 'Download and Configure Topbar' link – File download will be triggered.




Image 3: Topbar clinic view screen

Download and then run the download installer 'DesktopAppSetup1.x.x.x.exe'.
Self-extracting zip will run.




Image 4: Desktop App Installer

The installer will run




Image 5: Desktop App Setup

After the installation is complete, Topbar will start and the logon will be displayed on top of your screen.



Image 6: Topbar Logon