


Installation Preparation

1. Go to <https://topbar-users.pencs.com.au>
2. Sign in with the username and password provided to you or use your existing CAT account



The screenshot shows the TOPBAR Sign In page. At the top left is the TOPBAR logo. The page has a light gray header bar. The main content area is white. The title 'Sign In' is centered at the top of the main area. Below the title are two input fields: 'User Name' with a person icon and 'Password' with a lock icon. A blue 'LOGIN' button is centered below the fields. At the bottom left, there is a link 'Forgotten your password?' with a question mark icon.

TOPBAR

Sign In

User Name

Password

LOGIN

[? Forgotten your password?](#)

Once you have signed in to your account you will see your clinic listed and you can add other clinics if you work in more than one location. Clicking on the your clinic name on the tile will take you to the main Topbar configuration screen which can also be accessed at any by going to <https://topbar-users.pencs.com.au>

Welcome, Matthias.

You are currently a member of 2 clinics.



TEST CLINIC 2

1 New Street, Bowen Hills, 2482



MATTHIAS TEST CLINIC

1 New Street, Bowen Hills, 2482

Clicking on the your clinic name on the tile will take you to the Topbar clinic view screen

The clinic view screen below is described in more detail in the Topbar User guide.

MY CLINICS / MATTHIAS TEST CLINIC



VIEW CLINIC SETTINGS

Settings such as clinic hostname...



CONNECT TO CLINIC

Access your clinic




DOWNLOAD AND CONFIGURE TOPBAR

Install this on the users' PCs.
Version 2.1.22.3640

CLINIC SERVER
ALREADY INSTALLED

The Topbar Clinic Service is already installed on Win7VM and last contacted at 23/01/2018 4:05:00 PM
Version 2.1.22.3640



MATTHIAS TEST CLINIC

1 New Street,Bowen Hills,2482
Phone: 999999999

MANAGE USER LIST