De-identifying an Extract

1. Open CAT4 and select an extract that you would like to send to PAT CAT. Please ensure that you select an Extract which does NOT have the text "DEIDENTIFIED" or "DEIDENTIFIED_GTAGGED" in the extract name.



1. Open the Tools menu and select Deidentify Dataset



1. Press the OK button to confirm you wish to de-identify the extract



1. Press the OK button when the Success dialog appears



1. Notice that in the extracts panel, a new extract is listed with the text "DEIDENTIFIED" in the name. This is the file that is to be imported into PAT CAT.

PCS Clinical Audit Tool - Classic CAT			
File Edit View Tools Data Submission Pron	mpts Help		
Long Hot Colect Extension View Filter Report Population	The Case Case Case Case Case Case Case Cas	21/05/2014 4:25 PM	Char Films Restore
21/05/2014 4/29 PM 28 DEIDENTIFIED_MH	Demographics Data Quality Allergies Smoking	Alcohol Measures Pathology Disease Screening Medications Diabet	es SIP Items CKD Musculoskeletal CV Event Risk Immunisations Standard Reports MBS Items Maternal Health Care Visits Timeline
21/05/2014 4/29 PM 28 MHCHOPC001 21/05/2014 1/45 PM 29 DEIDENTIFIED GT	Select All Show Total Counts		Export Age bracket 5 😩 Print
21/05/2014 1:45 PM 28 MHCHOPC001			
21/05/2014 1:41 PM 28 DEIDENTIFIED_GT		Demographic Brea	kdown by Age [population = 28]
21/05/2014 1:41 PM 28 MHCHOPC001			
21/05/2014 12:29 PM 32 BP5000000		Females	= 12, Males = 16, Other = 0
1/04/2014 1:23 PM 3968 DEIDENTIFIED_SMC			
6/02/2014 12:13 PM 28 DEIDENTIFIED_MH		Earraia Maia	
1/01/2014 5:55 PM 20308 DEIDENTIFIED_8P	_	Penale Male	
	100+		
	95 - 99		
	90 - 94		
	85 - 89		
	80 - 84		
	75 - 79		
	/0 - /4		
	60-64-		
	2 55-59-		
	ž 50-54		
	45-49 —		
	40 - 44		
	35 - 39		
	30 - 34		
	25-29-		
	15-19-	· · · · · · · · · · · · · · · · · · ·	
	10 - 14		
	5-9-		
	0 - 4		
CLINICAL AUDIT TOOL		3 2 1	
		Nurrber of Pašerta	Number of Palients

More information on the De-identification and G-tagging process can be found in Appendix A.