Setting up CAT4 for Scheduler

The CAT4 scheduler preferences allow users to set the working folder. This folder is where CAT4 can update scheduled tasks and retrieve the scheduler history. It is monitored by the 'PEN CS Scheduler Service' so that scheduled tasks are run when required.

• CAT4 will automatically look for the default folder 'C:\Program Files\Pen Computer Systems\PEN CS Scheduler Service'. If this folder is found and the preferences working folder is empty CAT4 will pre-populate the preference. If you are working on a 64bit system the folder will be called C: \Program Files (x86)\Pen Computer Systems\PEN CS Scheduler Service

If you do not have Scheduler as a CAT4 menu option as shown below, you will need to configure the CAT4 Scheduler Preferences screen.



Follow the steps below to set up CAT4 for Scheduler

- 1. Open CAT4
- 2. Click the Edit > Preferences menu
- 3. From the Clinical Audit Preferences dialog choose the Scheduler tab
- 4. Note: If the Scheduler tab is not visible (due to your window size) a right arrow key will be available to scroll across. Alternatively, you can widen the window.

Scheduler Working Folder Directory should be automatically populated as C:*Program files(x86)\Pen Computer Systems\PCS Scheduler.* The specific **User Name** and **Password** for a practice is to be obtained from Pen CS *Support* for installation purpose only.

Clinical Audit Preferences	
Clinical Audit Best Practice Communicare Ferret Genie Mediflex Sidebar/Prompts Scheduler Data Submission T	Medical Director 3 Medinet MedTech32 MMEX practiX Zedmed
Scheduler Working Folder Directory	C:\Program Files (x86)\Pen Computer Systems\PCS Schedu
Scheduler Dedicated CAT Account Note	This is the account credentials that the Scheduled CAT collections will run under.
User Name	
Password	
	Validate Account
	OK Cancel

Once User Name and Password are filled in here, click the Validate Account button.

Clinical Audit Preferences		
Clinical Audit Best Practice Communicare Ferret Genie Mediflex Sidebar/Prompts Scheduler Data Submission T	Medical Director 3 Medinet MedTech32 MMEX practiX Zedmed	
Scheduler Working Folder Directory	C:\Program Files (x86)\Pen Computer Systems\PCS Schedu	
Scheduler Dedicated CAT Account Note	This is the account credentials that the Scheduled CAT collections will run under.	
User Name	scheduler.1000	
Password	Validate Account	
	OK Cancel	

Click OK to close the Scheduled Collection message box, which displays that the Account credentials verified and saved successfully.

If the account details can't be verified please contact Pen CS support.



Click OK for the Clinical Audit Preferences dialog box to close.

Clinical Audit Preferences	
Clinical Audit Best Practice Communicare Ferret Genie Mediflex Sidebar/Prompts Scheduler Data Submission	Medical Director 3 Medinet MedTech32 MMEX practiX Zedmed
Scheduler Working Folder Directory	C:\Program Files (x86)\Pen Computer Systems\PCS Schedu
Scheduler Dedicated CAT Account Note	This is the account credentials that the Scheduled CAT collections will run under.
User Name	scheduler.1000
Password	•••••
	Validate Account