

Choosing a preferred Data Extract Directory

If there are more than one CAT4 user in the practice then choose a **network drive** that is accessible to all. We recommend to use a network or shared drive that is included in your routine backups even if there is only one CAT4 user. If there are more than one user ensure that the selected location is accessible to all users.
Use the 'Edit > Preferences' option to set the same 'Extract Directory' for each user.

Clinical Audit Preferences

MedTech32MMEXpractiXStatHealthZedmedMediflexMedilinkThe PracticeSchedulerData SubmissionTopbar

Clinical AuditRefugee HealthBest PracticeCommunicareFerretGenieGPCompleteMedical Director 3Medinet

GP Application

Medical Director 3

Extract Directory

C:\ClinicalAudit

...

Saved Filters Directory

C:\ClinicalAudit\SavedFilters

...

Saved Program Views Directory

C:\ClinicalAudit\SavedProgramViews

...

Saved Views Directory

C:\ClinicalAudit\SavedViews

...

Practice ID

Practice Name

Practice Postcode

Patient Consent Withdrawn File Directory

C:\ClinicalAudit\Settings

...

Patient Consent Withdrawn File Directory Note

It is recommended that this directory is set to a shared folder that all CAT users have access to.

Run CAT in 64 bit mode

Yes

OK

Cancel



More than 1 user can work with the same snapshot at the same time.

All directories should be configured to be on the server where CAT4 is installed and extracts are saved. Please make sure that all CAT4 users have access to the extract folder, as otherwise no data will be visible on the workstation.

Genie/practiX users will need to make sure this folder matches the folder they select from the Genie/practiX application when they do the extract.