

Medical Director Requirements

Select your Database Configuration from the drop-down list. These will match the configurations you have available through the 'HCN Maintenance' icon on your desktop. You can also change the settings for extracting archived and deceased patients which by default is set to 'No'. This is generally not needed in General Practice but is required for the [NKPI report](#). The Location list settings is relevant for clinics that share one server for multiple locations - more details are described here: [CAT4 MD3 Collection by Location](#)

Please note that the "Progress Note Batch Size" settings are only required if there is a problem with corrupt tables in your MD3 installation. Please contact Pen CS support to assist if you are experiencing issues with no active patients showing in your CAT4.

Clinical Audit Preferences

MedTech32 | MMEX | practiX | Zedmed | Mediflex | Medilink | The Practice | Scheduler | Data Submission | Topbar
Clinical Audit | Refugee Health | Best Practice | Communicare | Ferret | Genie | GPComplete | **Medical Director 3** | Medinet

Local AZdex Directory: C:\Program Files (x86)\Health Communication Network\Med ...

Database Configuration: MD Live Data

Database Server (optional override):

Database Server Note: If you change the server you will need to restart CAT to refresh the Location List.

Document Server: C:\ManageEzy\

Document Library: HCN

Billing Software: Pracsoft 3

Mediflex/Medilink/The Practice Note: If you are using Medical Director 3 and Mediflex, Medilink or The Practice, please check the relevant tab to ensure the 'Application Directory' Setting is correct.

Progress Note Batch Size (optional override):

Progress Note Batch Size Note: This setting is only required where the extract contains no Active patients. This is caused by a Progress Notes table indexing error that can be resolved using batching.

Extract Archived and Deceased Patients: No

Practice Location List (comma separated)
(Enter ALL for a single extract, IND for an extract per Pracsoft location)

Program Location List (comma separated)
Prefix entered in patient record no
(Enter ALL for a single extract)

OK Cancel