


# Spirometry

Measurements > Spirometry displays your data as a breakdown of the ratio of POST FEV1/FVC.  
FEV1 is the forced expired volume in one second  
FVC is the forced vital capacity  
FEV1/FVC provides an index of airflow limitation; a value < 0.7 indicates limitation  
The breakdown is shown as

- POST FEV1/FVC < 0.07
- POST FEV1/FVC >= 0.07
- Not recorded



## Spirometry Mappings with clinical systems

The Spirometry mappings are provided in the 'General Data Mappings' Appendices under Measurements for each clinical application.